## ***Anti- GluA2 (GluR2)***

*(AMPA-type glutamate receptor subunit-2)*

**Code Number** : GluR2C-Rb-Af1050 (rabbit, RRID : AB\_2571754)

**Size** : 20 µg and 50 µg / See label on vial  
(affinity-purified with antigen polypeptide)


**Formulation** : Liquid ; 200 µg/ml in PBS with 0.05% NaN<sub>3</sub>.

**Storage** : Store at cool temp. (2-10°C)

The antibody can be stored at 2-10°C. The antibody can be also aliquotted and stored at -80°C for long-term storage. Avoid repeated freeze-thawing. Non-hazardrous. No MSDS required.

**Species** : rabbit, polyclonal

**Antigen** : mouse GluR2, C-terminal  
847-863 aa  
(X57498 ; KVAKNAQNINPSSSQNS)

**Specificity** : mouse (others not tested)

Immunoblot detects a single protein band at 100kDa, with no cross reactivity to other iGluR subunits, including GluR3.

**Applications** : In general, affinity-purified antibody is used at around 1 microgram/ml for immunoblot and immunohistochemistry. The most appropriate concentration should be determined by users, because it depends on contents in given cells, tissues and organs.

**Research Use** : For research use only, not for use in diagnostic procedures.

**Remarks** : For immunohistochemistry for neuronal iGluRs, users should adopt postembedding immunogold for electron microscopic detection and protease predigestion for light microscopic detection (see the below reference). For glial GluR, these antigen-exposing methods are not necessary (unpublished information).

**Reference** : 1) Shimuta, M., Yoshikawa, M., Fukaya, M., Watanabe, M., Takeshima, H., Manabe, T. (2001) Postsynaptic modulation of AMPA receptor-mediated synaptic responses and LTP by

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Distributor : NITTOBO MEDICAL CO., LTD
Manufacture : FRONTIER INSTITUTE CO., LTD