

## *Anti-Acetyl-Spermine mAb*

### Product Information:

---

1) <b>Catalog #:</b>	MSFR107290
<b>Item Name:</b>	Acetyl-Spermine mAb 75-20
<b>Size :</b>	50µg (Protein G purified)
<b>Species :</b>	Mouse
<b>Product Code :</b>	A-SPERMINE-MO-75-20
<b>RRID :</b>	AB_2848124

---

**Formulation:** PBS solution, pH7.4

A-SPERMINE-MO-75-20 : 0.14 mg/mL


**Storage:** Store at -20 °C. The antibody can be also aliquoted and stored at -80 °C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous.

**Species:** mouse, monoclonal

**Antigen:** Acetyl-Spermine

**Reactivity:**

**Indirect competitive STD curve /N- AcSPM**



**Anti-Acetyl-Spermine mAb cross reactivity \***

Ac-SPM	100%
N8-AcSPD	<0.1%
N <sub>1</sub> N <sub>12</sub> Di Ac SPM	230%
N <sub>1</sub> N <sub>8</sub> Di Ac SPD	0.90%

\* Use Ac SPM conjugate coated plate

**Applications:** ELISA

**Research Use:** For research use only