

Anti-CD224 antibody

Purified antibody

(CD224, GGT1 gamma-glutamyltransferase 1)

Product Information:

1) Catalog #:	MSFR106750
Item Name:	CD224 mAb FRT39 (purified) 100ug
Size :	100µg (affinity-purified, protein G)
Species :	Mouse
Product Code :	CD224-Mo-MY035
RRID :	AB_2750921

Clone name : FRT39

Other name : GGT, GTG, GGTD, CD224, GGT 1, D22S672, D22S732

Size : 100 µg (purified , proteinG)

Formulation : liquid. (1.0 mg / ml in PBS)

Storage : The antibody are refrigerated (2-10°C) at the time of shipment, but should be divided into small portions and frozen (-80 °C) after arrival. This antibody can stored at -80°C for long-term storage. Avoid repeated freeze-thaw cycles. Non-hazardous.


Species : mouse, monoclonal

Isotype : IgG1, κ


Antigen : MIA-PaCa2/Human

Specificity : human (others not tested)

This antibody is specific to CD224.



cell line:HCC827
2nd antibody:anti-mouse IgG FITC



cell line:HS-Sultan
2nd antibody:anti-mouse IgG FITC

Applications : Flow Cytometry (FC), ImmunoPrecipitation (IP)

Reseaech Use : For research use only, not for use in diagnostic procedures.

Remark : FC recommended concentration

(cell conc. 5×10^5 cell / test & antibody conc. $0.5 \mu\text{g}$ / test)

Reference :