

Anti-CD224 antibody

Purified antibody

(CD224, GGT1 gamma-glutamyltransferase 1)

code Number : CD224-Mo-MY035 (RRID : AB_2750921)

Clone name : FRT39

Other name : GGT, GTG, GGTD, CD224, GGT 1, D22S672, D22S732

Size : 100 µg (purified , proteinG)

Formulation : liquid. (1.10 mg / ml in PBS with 0.03% ProClin.)

Storage : Store at cool temp (2-10°C).

The antibody can be stored at 2-10°C. This antibody can also be aliquotted and stored at -80°C for long-term storage. Avoid repeated freeze-thaw cycles. Non-hazardous. No MSDS required.


Species : mouse, monoclonal

Isotype : IgG1, κ


Antigen : MIA-PaCa2/Human

Specificity : human (others not tested)

This antibody is specific to CD224.



cell line:HCC827
2nd antibody:anti-mouse IgG FITC



cell line:HS-Sultan
2nd antibody:anti-mouse IgG FITC

ご注意 : 本商品には 0.03 % の ProClin が入っています。 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD
TEL : +81-(0)3-4582-5451
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD

Applications : Flow Cytometry (FC), ImmunoPrecipitation (IP)

Research Use : For research use only, not for use in diagnostic procedures.

Remark : FC recommended concentration

(cell conc. 5×10^5 cell / test & antibody conc. 0.5 μg / test)

Reference :

ご注意：本商品には 0.03 % の ProClin が入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD
TEL : +81-(0)3-4582-5451
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD