

Anti-CD63 Monoclonal Antibody

Purified antibody

(LAMP3, Lysosomal-Associated Membrane Protein 3)

Product Information:

1) Catalog #:	MSFR100750
Item Name:	CD63 mAb FRT8 (purified) 100ug
Size :	100µg (affinity-purified, protein G)
Species :	Mouse
Product Code :	CD63-Mo-MY006
RRID :	AB_2722743
2) Catalog #:	MSFR107340
Item Name:	CD63 mAb FRT8 (purified) 1 mg
Size :	1mg (affinity-purified, protein G)
Species :	Mouse
Product Code :	CD63-Mo-MY006
RRID :	AB_2722743

Clone name : FRT8

Other name : LAMP3, Lysosomal-Associated Membrane Protein 3

Formulation : liquid. (1 mg / ml in PBS)

Storage : Store at 4°C. This antibody can also be aliquoted and stored at -80°C for long-term storage. Avoid repeated freeze-thaw cycles. Non-hazardous. No MSDS required.


Species : mouse, monoclonal

Isotype : IgG2a, κ


Antigen : MIA-PaCa2/Human

Specificity : human (others not tested)

This antibody is specific to CD63.



cell line:PC9
2nd antibody:anti-mouse IgG FITC



cell line:A375
2nd antibody:anti-mouse IgG FITC

Applications : Flow Cytometry (FC), Immunoprecipitation (IP)

Research Use : For research use only, not for use in diagnostic procedures.

Remark : FC recommended concentration

(cell conc. 5×10^5 cell / test & antibody conc. $0.5 \mu\text{g}$ / test)

Reference :