

Anti-CD9 Monoclonal Antibody

Purified antibody

(TSPAN-29)

Product Information:

1) Catalog #:	MSFR106710
Item Name:	CD9 mAb FRT3 (purified) 100ug
Size :	100µg (affinity-purified, protein G)
Species :	Mouse
Product Code :	CD9-Mo-MY050
RRID :	AB_2754556
2) Catalog #:	MSFR107360
Item Name:	CD9 mAb FRT3 (purified) 1 mg
Size :	1mg (affinity-purified, protein G)
Species :	Mouse
Product Code :	CD9-Mo-MY050
RRID :	AB_2754556

Clone name : FRT3

Other name : MIC3; MRP-1; BTCC-1; DRAP-27; TSPAN29; TSPAN-29

Formulation : liquid. (1 mg / ml in PBS)

Storage : The antibody are refrigerated (2-10°C) at the time of shipment, but should be divided into small portions and frozen (-80 °C) after arrival. This antibody can store at -80°C for long-term storage. Avoid repeated freeze-thaw cycles. Non-hazardous.

Species : mouse, monoclonal

Isotype : IgG2b, κ

Antigen : MIA-PaCa2/Human


ご注意：本商品には 0.03 % の ProClin が入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD
TEL : +81-(0)3-4582-5451
email : nmd-tokuhan@nittobogrp.com


Manufacture : FRONTIER INSTITUTE CO., LTD

Specificity : human (others not tested)

This antibody is specific to CD9.



cell line:PC9
2nd antibody:anti-mouse IgG FITC



cell line:A375
2nd antibody:anti-mouse IgG FITC

Applications : Flow Cytometry (FC), Immunoprecipitation (IP)

Research Use : For research use only, not for use in diagnostic procedures.

Remark : FC recommended concentration

(cell conc. 5×10^5 cell / test & antibody conc. $0.5 \mu\text{g}$ / test)

Reference :

ご注意 : 本商品には 0.03 % の ProClin が入っています。 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD
TEL : +81-(0)3-4582-5451
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD