

Anti-DsRed (red fluorescent protein)

Product Information:

1)	Catalog #:	MSFR101390
	Item Name:	DsRED pAb (Rb) 20ug
	Size :	20µg (affinity-purified with antigen polypeptide)
	Species :	Rabbit
	Product Code :	DsRED-Rb-Af710
	RRID :	AB_2571647
2)	Catalog #:	MSFR101400
	Item Name:	DsRED pAb (Rb) 50ug
	Size :	50µg (affinity-purified with antigen polypeptide)
	Species :	Rabbit
	Product Code :	DsRED-Rb-Af710
	RRID :	AB_2571647
3)	Catalog #:	MSFR101410
	Item Name:	DsRED pAb (GP) 20ug
	Size :	20µg (affinity-purified with antigen polypeptide)
	Species :	Guinea pig
	Product Code :	DsRED-GP-Af360
	RRID :	AB_2571648
4)	Catalog #:	MSFR101420
	Item Name:	DsRED pAb (GP) 50ug
	Size :	50µg
	Species :	Guinea pig
	Product Code :	DsRED-GP-Af360
	RRID :	AB_2571648

Formulation : Liquid ; 200 µg/ml in PBS with 0.05% NaN₃.
(affinity-purified with antigen polypeptide)

Storage : Store at cool temp. (2-10°C)

The antibody can be stored at 2-10°C. The antibody can be also aliquoted and stored at -80°C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous. No MSDS required.

Species : rabbit / guinea pig, polyclonal

Antigen : DsRED, 1-225 aa (GU253312)

Specificity : mouse (others not tested)

This selectively stains cells in DsRED-transgenic mice, but not wild-type mice.

Applications : In general, affinity-purified antibody is used at around 1 microgram/ml for immunoblot and immunohistochemistry. The most appropriate concentration should be determined by users, because it depends on contents in given cells, tissues and organs.

ご注意 : 本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Research Use : For research use only, not for use in diagnostic procedures.

Remarks :

Reference : paper in submission

ご注意：本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD
TEL : +81-(0)3-4582-5451
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD