

## *Anti-EpCAM antibody*

### *Purified antibody*

*(Epithelial Cell Adhesion Molecule, CD326)*

**code Number** : EpCAM-Mo-MY012 (RRID : AB\_2750573)

**Clone name** : FRT16

**Other name** : Epithelial Cell Adhesion Molecule, CD326

**Size** : 100 µg (purified , proteinG)

**Formulation** : liquid. ( 1 mg / ml in PBS with 0.03% ProClin.)

**Storage** : Store at cool temp (2-10°C).

The antibody can be stored at 2-10°C. This antibody can also be aliquotted and stored at -80°C for long-term storage. Avoid repeated freeze-thaw cycles. Non-hazardous. No MSDS required.


**Species** : mouse, monoclonal

**Isotype** : IgG1, κ


**Antigen** : NCI-H322/Human

**Specificity** : human (others not tested)

This antibody is specific to EpCAM.



cell line:PC9  
2<sup>nd</sup> antibody:anti-mouse IgG FITC



cell line:A375  
2<sup>nd</sup> antibody:anti-mouse IgG FITC

ご注意 : 本商品には 0.03 % の ProClin が入っています。 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD  
TEL : +81-(0)3-4582-5451  
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD

**Applications :** Flow Cytometry (FC), ImmunoPrecipitation (IP)

**Research Use :** For research use only, not for use in diagnostic procedures.

**Remark :** FC recommended concentration

(cell conc.  $5 \times 10^5$  cell / test & antibody conc. 0.5  $\mu$ g / test)

**Reference :**

ご注意：本商品には 0.03 % の ProClin が入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD  
TEL : +81-(0)3-4582-5451  
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD