

Anti-GAT3

(Plasmalemmal GABA transporter-3, Slc6a11)

Product Information:

1) Catalog #:	MSFR106530
Item Name:	GAT3(1-30) pAb (Rb) 20ug
Size :	20µg (affinity-purified with antigen polypeptide)
Species :	Rabbit
Product Code :	GAT3(1-30)-Rb-Af630
RRID :	AB_2728722
2) Catalog #:	MSFR106540
Item Name:	GAT3(1-30) pAb (Rb) 50ug
Size :	50µg (affinity-purified with antigen polypeptide)
Species :	Rabbit
Product Code :	GAT3(1-30)-Rb-Af630
RRID :	AB_2728722
3) Catalog #:	MSFR106550
Item Name:	GAT3(590-627) pAb (Rb) 20ug
Size :	20µg (affinity-purified with antigen polypeptide)
Species :	Rabbit
Product Code :	GAT3(590-627)-Rb-Af440
RRID :	AB_2728731
4) Catalog #:	MSFR106560
Item Name:	GAT3(590-627) pAb (Rb) 50ug
Size :	50µg (affinity-purified with antigen polypeptide)
Species :	Rabbit
Product Code :	GAT3(590-627)-Rb-Af440
RRID :	AB_2728731

Formulation : Liquid ; 200 µg/ml in PBS with 0.05% NaN₃.

Storage : Store at cool temp. (2-10°C)

The antibody can be stored at 2-10°C. The antibody can be also aliquoted and stored at -80°C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous. No MSDS required.

Species : rabbit, polyclonal


Antigen : mouse, 1-30 aa, 590-627aa (NM_172890.3)

ご注意 : 本商品には 0.1 % 未満のアジ化ナトリウムが入っています。 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Specificity : mouse (others not tested)

Immunoblot detects a protein band at 60 kDa in mouse brains.

Applications : This can be used for immunohistochemistry and immunoblot.



Research Use : For research use only, not for use in diagnostic procedures.

Remarks : Note that GAT3 (Slc6a11) is named GAT4 in mice. Mouse GAT3 corresponds to Slc6a13.

Reference :

ご注意 : 本商品には 0.1 % 未満のアジ化ナトリウムが入っています。 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD
TEL : +81-(0)3-4582-5451
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD