

## *Anti-GAT-3*

*(Plasmalemmal GABA transporter-3, Slc6a11)*

**Code Number** : GAT3-Rb-Af630 (1-30) (rabbit, RRID : AB\_2728722)  
 : GAT3-Rb-Af440 (590-627) (rabbit, RRID : AB\_2728731)

**Size** : 20 µg and 50 µg / See label on vial  
 (affinity-purified with antigen polypeptide)

**Formulation** : Liquid ; 200 µg/ml in PBS with 0.05% NaN<sub>3</sub>.

**Storage** : Store at cool temp. (2-10°C)

The antibody can be stored at 2-10°C. The antibody can be also aliquotted and stored at -80°C for long-term storage. Avoid repeated freeze-thawing. Non-hazardrous. No MSDS required.


**Species** : rabbit, polyclonal

**Antigen** : mouse, 1-30 aa, 590-627aa (NM\_172890.3)

**Specificity** : mouse (others not tested)

Immunoblot detects a protein band at 60 kDa in mouse brains.

**Applications** : This can be used for immunohistochemistry and immunoblot.



**Research Use** : For research use only, not for use in diagnostic procedures.

ご注意 : 本商品には 0.1 % 未満のアジ化ナトリウムが入っています。 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD  
 TEL : +81-(0)3-4582-5451  
 email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD

**Remarks** : Note that GAT3 (Slc6a11) is named GAT4 in mice. Mouse GAT3 corresponds to Slc6a13.

**Reference** :

ご注意：本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD  
TEL : +81-(0)3-4582-5451  
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD