

Anti-GST tag

Product Information:

1) Catalog #:	MSFR080670
Item Name:	Anti GST-Tag antibody
Size :	100µg (Affinity purified)
Species :	Rabbit
Product Code :	GSTAB


Formulation: PBS solution with 0.03% ProClin as preservative.

Storage: Store at 4 °C. Avoid repeated freeze-thawing. Non-hazardous.

Species: rabbit, polyclonal

Antigen: GST tagged recombinant protein

Specificity: The antibody reacts with GST fusion protein



Applications: Western blotting (x1/500 – x1/3000, Optimal dilution/concentration should be determined by the researcher.)

Not tested in other applications.

Research Use: For research use only, not for use in diagnostic procedures.