

Anti-Human IgG4 mAb

Product Information:

1) Catalog #:	MSFR107270
Item Name:	Human IgG4 mAb 16-3
Size :	50µg (Protein G purified)
Species :	Mouse
Product Code :	HIGG4-MO-16-3
RRID :	AB_2848122
2) Catalog #:	MSFR107280
Item Name:	Human IgG4 mAb 19-3
Size :	50µg (Protein G purified)
Species :	Mouse
Product Code :	HIGG4-MO-19-3
RRID :	AB_2848123

Formulation: PBS solution, pH7.4

- 1) HIGG4-MO-16-3: 0.2 mg/mL
- 2) HIGG4-MO-19-3: 0.09 mg/mL

Storage: Store at -20 °C. The antibody can be also aliquoted and stored at -80 °C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous.

Species: mouse, monoclonal IgG4


Antigen: Human IgG4

ご注意：本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。


Distributor : NITTOBO MEDICAL CO., LTD Manufacture : FRONTIER INSTITUTE CO., LTD
TEL : +81-(0)3-4582-5451 email : nmd-tokuhan@nittobogrp.com

Cross reactivity: The unit of concentration: (ng/mL).

Anti Human IgG4, #16-3



Anti Human IgG4, #19-3



Applications: ELISA, WB

Research Use: For research use only

ご注意：本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。