

Anti-human β 2 microglobulin mAb

Product Information:

1) Catalog #:	MSFR107310
Item Name:	Human Beta-2 microglobulin mAb 13-5
Size :	50 μ g ProteinG purified
Species :	Mouse
Product Code :	B2MG-MO-13-5
RRID :	AB_2858224
2) Catalog #:	MSFR107320
Item Name:	Human Beta-2 microglobulin mAb 21-3
Size :	50 μ g ProteinG purified
Species :	Mouse
Product Code :	B2MG-MO-21-3
RRID :	AB_2858225
3) Catalog #:	MSFR107330
Item Name:	Human Beta-2 microglobulin mAb 20-3
Size :	50 μ g ProteinG purified
Species :	Mouse
Product Code :	B2MG-MO-20-3
RRID :	AB_2858226

Formulation: PBS solution, pH7.4


1) B2MG-MO-13-5 :	1 mg/mL
2) B2MG-MO-21-3 :	1 mg/mL
3) B2MG-MO-20-3 :	0.4 mg/mL

Storage: Store at -20 °C. The antibody can be also aliquoted and stored at -80 °C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous.

Species: mouse, monoclonal

Antigen: Human β 2 microglobulin

Cross reactivity: Human serum albumin: < 0.1 %
Human α 2 macroglobulin: < 0.1 %
Human γ globulin: < 0.1 %



Applications: ELISA, WB, IP

Research Use: For research use only