

## ***Anti-Kv1.2***

*(Shaker-type potassium channel, type 2)*

**Code Number** : Kv1.2-Rb-Af240 (rabbit, RRID : AB\_2571789)

**Size** : 20 µg and 50 µg / See label on vial  
(affinity-purified with antigen polypeptide)

**Formulation** : Liquid ; 200 µg/ml in PBS with 0.05% NaN<sub>3</sub>.

**Storage** : Store at cool temp. (2-10°C)

The antibody can be stored at 2-10°C. The antibody can be also aliquotted and stored at -80°C for long-term storage. Avoid repeated freeze-thawing. Non-hazardrous. No MSDS required.


**Species** : rabbit , polyclonal

**Antigen** : mouse, 482-496 aa (NM\_008417)

**Specificity** : mouse (others not tested)

Immunoblot detects a protein band  
at 90 kDa In mouse, brains.

**Applications** : This can be used for  
immunohistochemistry and immunoblot.



**Research Use** : For research use only, not for use in diagnostic procedures.

**Remarks** :

**Reference** : Iwakura A, Uchigashima M, Miyazaki T, Yamasaki M, Watanabe M: Lack of molecular-anatomical evidence for GABAergic influence upon axon initial segment of cerebellar Purkinje cells by the pinceau formation. **J. Neurosci.** 32:9438-9448.

ご注意 : 本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD  
TEL : +81-(0)3-4582-5451  
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD