

## *Anti-NCBE*

*(Na<sup>+</sup>-driven Cl<sup>-</sup>/HCO<sub>3</sub><sup>-</sup> exchanger pr Slc4a10)*

### Product Information:

---

1) Catalog #:	MSFR104300
Item Name:	NCBE (Slc4a10) pAb (Rb) 20ug
Size :	20µg (affinity-purified with antigen polypeptide)
Species :	Rabbit
Product Code :	NCBE-Rb-Af760
RRID :	AB_2571807
2) Catalog #:	MSFR104310
Item Name:	NCBE (Slc4a10) pAb (Rb) 50ug
Size :	50µg (affinity-purified with antigen polypeptide)
Species :	Rabbit
Product Code :	NCBE-Rb-Af760
RRID :	AB_2571807
3) Catalog #:	MSFR104320
Item Name:	NCBE (Slc4a10) pAb (GP) 20ug
Size :	20µg (affinity-purified with antigen polypeptide)
Species :	Guinea pig
Product Code :	NCBE-GP-Af870
RRID :	AB_2571808
4) Catalog #:	MSFR104330
Item Name:	NCBE (Slc4a10) pAb (GP) 50ug
Size :	50µg (affinity-purified with antigen polypeptide)
Species :	Guinea pig
Product Code :	NCBE-GP-Af870
RRID :	AB_2571808

---

**Formulation** : Liquid ; 200 µg/ml in PBS with 0.05% NaN<sub>3</sub>.  
(affinity-purified with antigen polypeptide)

**Storage** : Store at cool temp (2-10°C).

The antibody can be stored at 2-10°C. The antibody can be also aliquoted and stored at -80°C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous. No MSDS required.

**Species** : rabbit / guinea pig, polyclonal

**Antigen** : mouse NCBE, 17-47 aa (AB033759)

**Specificity** : mouse (others not tested)

Immunoblot detects a single protein band at 130-140kDa in brain homogenates and HEK293T cells expressing Slc4a10. Immunocytochemistry selectively labels HEK293T cells expressing Slc4a10, but not untransfected cells. Preabsorption test abolishes immunohistochemical labeling in adult mouse brains, as reported in Song et al. (2014).

ご注意 : 本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

**Applications** : In general, affinity-purified antibody is used at around 1 microgram/ml for immunoblot and immunohistochemistry. The most appropriate concentration should be determined by users, because it depends on contents in given cells, tissues and organs.

**Research Use** : For research use only, not for use in diagnostic procedures.

**Remarks** : Rabbit and guinea pig antibodies are similar in titer and specificity.

**Reference** : Song X, Yamasaki M, Miyazaki T, Konno K, Uchigashima M, Watanabe M: Neuron type- and input pathway-dependent expression of Slc4a10 in adult mouse brains. **Eur J Neurosci**, 40:2797-2810, 2014.

ご注意：本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD  
TEL : +81-(0)3-4582-5451  
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD