

Anti-RyR1

(ryanodine receptor 1, skeletal type)

Product Information:

1) Catalog #:	MSFR107140
Item Name:	RyR1 pAb (Rb)
Size :	20µg (affinity-purified with antigen polypeptide)
Species :	Rabbit
Product Code :	RyR1-Rb-Af660
RRID :	AB_2858216
2) Catalog #:	MSFR107150
Item Name:	RyR1 pAb (Rb)
Size :	50µg (affinity-purified with antigen polypeptide)
Species :	Rabbit
Product Code :	RyR1-Rb-Af660
RRID :	AB_2858216

Formulation: Liquid; 200 µg/ml in PBS with 0.05% NaN₃.

Storage: Store at cool temperature (2-10°C).

The antibody can be stored at 2-10°C. The antibody can be also aliquoted and stored at -80°C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous. No MSDS required.

Species: rabbit, polyclonal

Antigen: mouse, 1343-1404 aa (AY268935)

Specificity: mouse (others not tested)

This antibody stains skeletal, but not cardiac, muscle. In the brain, immunohistochemical labeling is particularly intense in the cerebellar cortex, where aggregated puncta are crowded in Purkinje cell perikarya and dendrites.

ご注意：本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

