

Anti-Tubulin β -III

Code Number : TubIII-Rb-Af220 (rabbit, RRID : AB_2571848)
: TubIII-GP-Af760 (guinea pig, RRID : AB_2571849)

Size : 20 μ g and 50 μ g / See label on vial
(affinity-purified with antigen polypeptide)

Formulation : Liquid ; 200 μ g/ml in PBS with 0.05% NaN₃.


Storage : Store at cool temp (2-10°C).

The antibody can be stored at 2-10°C. The antibody can be also aliquotted and stored at -80°C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous. No MSDS required.

Species : rabbit / guinea pig, polyclonal

Antigen : mouse TubIII, C-terminal 17 aa
(BC088749)

Specificity : mouse (others not tested)
This selectively stains neuronal dendrites,
perikarya, and axons.



Applications : In general, affinity-purified antibody is used at around 1 microgram/ml for immunoblot and immunohistochemistry. The most appropriate concentration should be determined by users, because it depends on contents in given cells, tissues and organs.

Research Use : For research use only, not for use in diagnostic procedures.

Remarks : Rabbit and guinea pig antibodies are similar in the titer and specificity.

Reference :

ご注意 : 本商品には 0.1 % 未満のアジ化ナトリウムが入っています。誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD
TEL : +81-(0)3-4582-5451
email : nmd-tokuhan@nittobogrp.com

Manufacture : FRONTIER INSTITUTE CO., LTD