

Goat IgG , purified


Product Information:

1) Catalog #:	MSFR070020
Item Name:	Goat IgG purified 500 μ g
Size :	500 μ g (Protein G purified)
Species :	Goat
Product Code :	7002
RRID :	-

Formulation: Liquid; 10mg/mL in PBS (antiseptic free)

Storage: Store at -20 °C. The antibody can be also aliquoted and stored at -80 °C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous. No MSDS required.

Species: goat, polyclonal



12% SDS PAGE , Sample 5 μ g

Applications: Blocking, Antigen, Negative control

Research Use: For research use only, not for use in diagnostic procedures.