

## *Anti - Plasmodium LDH (pLDH) , polyAb*

**Code Number** : PLDH-CH-81

**Size** : 100µg or order in bulk

**Formulation** : PBS solution , Affinity purified


**Storage** : Store at -20 °C. Avoid repeated freeze-thawing. Non-hazardrous. No MSDS required.

**Species** :Chicken, polyclonal, IgY ( Egg yolk )

**Antigen** :

*Plasmodium rLDH*

**Specificity**: *p LDH, N terminal*



*Western blot of P.falciparum rLDH*

**Applications** : User optimized (ELISA, Lateral Flow, Western blotting)

**Research Use** : For research use only.

**Remarks** :

ご注意 : 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD

Manufacture : FRONTIER INSTITUTE CO., LTD

TEL : +81-(0)3-4582-5451

email : nmd-tokuhan@nittobogrp.com