

Anti – Plasmodium LDH (pLDH), mAb

Code Number : PLDH-MO-704

Size : 50µg or order in bulk

Formulation : PBS solution, ProteinG purified


Storage : Store at -20°C. Avoid repeated freeze-thawing. Non-hazardrous. No MSDS required.

Species : mouse, monoclonal, IgG1

Antigen :

Plasmodium LDH

Specificity : *pLDH, N terminal*



Western blot of P. berghei ANKA infected RBC lysate

Applications : ELISA, Western blotting, IC

Research Use : For research use only.

Remarks :

ご注意 : 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD

Manufacture : FRONTIER INSTITUTE CO., LTD

TEL : +81-(0)3-4582-5451

email : nmd-tokuhan@nittobogrp.com