

Rabbit IgG , purified


Product Information:

1) Catalog #:	MSFR070010
Item Name:	Rabbit IgG purified 500 μ g
Size :	500 μ g (ProteinG purified)
Species :	Rabbit
Product Code :	7001
RRID :	

Formulation: Liquid; 4mg/mL in PBS (antiseptic free)

Storage: Store at -20 °C. The antibody can be also aliquoted and stored at -80 °C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous. No MSDS required.

Species: rabbit, polyclonal



12% SDS PAGE , Sample 5 μ g

Applications: Blocking, Antigen, Negative control

Research Use: For research use only, not for use in diagnostic procedures.