

Rabbit IgG, purified


Code Number : 7001

Size : 500 μ g (ProteinG purified)

Formulation : Liquid; 4mg/mL in PBS (antiseptic free)

Storage : Store at -20 °C. The antibody can be also aliquotted and stored at -80 °C for long-term storage. Avoid repeated freeze-thawing. Non-hazardrous. No MSDS required.

Species : rabbit, polyclonal



Applications : Blocking, Antigen, Negative control

Research Use : For research use only, not for use in diagnostic procedures.

ご注意 : 誤って目や口に入ったり、皮膚に付着した場合は大量の水で洗い流してください。

Distributor : NITTOBO MEDICAL CO., LTD

Manufacture : FRONTIER INSTITUTE CO., LTD

TEL : +81-(0)3-4582-5451

email : nmd-tokuhan@nittobogrp.com